

MODULE 2: Depositing cells into the paper scaffolds




Figure 1. Demonstration of the deposition of cell-laden collagen into a A) Large Single Zone (3.14 mm^3 culture area) and B) Small Single Zone (0.22 mm^3 culture area). The scale bars represent A) 10.00 mm and B) 2.65 mm.

Materials:

- 1X DPBS
- 10 cm untreated petri dish
- 96-well plate
- High glucose DMEM with supplements
- T150 flask of fluorescently labeled HepG2 cells (transfected or labeled)
- TrypLE express

- Pre-sterilized small single zones
- 15 mL conical tube
- 1.5 mL microcentrifuge tube (Fisher Scientific, no. 05-408-129)
- 96-well plate

- Aven 3rd hand with tweezers (Amazon, no. B006RB93IG)
- Benchtop 15 mL centrifuge
- Biosafety cabinet
- Inverted microscope
- Cell counter
- Flatbed fluorescence scanner
- Ice bucket
- Inverted fluorescence microscope with camera
- Microcentrifuge

Micropipettes and pipette tips
pH paper, pH 6.0-8.0 (Fisher Scientific, no. 13-640-510)
Plate reader capable of fluorescence and chemiluminescence readouts
Pre-sterilized small paper scaffolds

Steps:

1. In the Biosafety cabinet, detach the fluorescently labeled HepG2 cells from the plate using TrypLe express and count following standard culture techniques. Record the cell concentration as cells/ mL (C_{stock}).

Note. If the cells don't constitutively express a fluorescent protein, load the cells with dye prior to depositing them into the paper scaffolds. Some stains and dyes will nonspecifically absorb to the paper fibers making analysis challenging.

2. Calculate the number of cells ($N_{total\ cells}$) needed to prepare $N_{scaffolds}$ containing 4.0×10^4 cells/scaffold with Eqn. 1.

$$N_{total\ cells} = N_{scaffolds} \times \frac{cells}{scaffold} \times 1.5 \quad \text{Eqn. 1}$$

Note. A factor of 1.5 is included to ensure there is enough of the cell hydrogel suspension when depositing cells into the paper scaffolds. The collagen suspension is viscous, pipette slowly to ensure transfer of the appropriate volume and avoid adding bubbles.

3. Determine the volume of the stock cell suspension required ($V_{Cell\ stock}$) required to obtain $N_{total\ cells}$ with Eqn. 2.

$$V_{Cell\ stock} = \frac{N_{total\ cells}}{\left(C_{stock} \frac{cells}{mL}\right)} \quad \text{Eqn. 2}$$

4. Calculate the volume of collagen needed to prepare both the cell suspension ($V_{Cell\ dep}$) and blank scaffolds ($V_{dep, blank}$) with Eqns. 3 and 4. The 1.5 factor is to match Eqn. 1 and ensure there is enough collagen to deposit cells into all the scaffolds needed for the experiment.

$$V_{Cell\ dep} = 1.5 \times (N_{scaffolds} \times 0.5\ \mu L) \quad \text{Eqn. 3}$$

$$V_{dep, blank} = 1.5 \times (N_{scaffolds} \times 0.5\ \mu L) \quad \text{Eqn. 4}$$

Note. $V_{dep, blank}$ are scaffolds that only contain collagen. These blank scaffolds are needed during analysis.

5. Calculate the volume of each component needed to prepare a 1.2 mg/mL solution of collagen I with Eqn. 5 or 6.

$$V_{total} = V_{Cell\ dep} + V_{dep, blank} \quad \text{Eqn. 5}$$

or

$$V_{total} = V_{stock\ collagen} + V_{PBS} + V_{NaOH} + V_{water} \quad \text{Eqn. 6}$$

In these equations, $V_{stock\ collagen}$ is the volume of collagen stock (Eqn. 7), V_{PBS} is the volume of a 10X PBS solution (Eqn. 8), V_{NaOH} is the volume of a 1 M solution of sodium hydroxide (Eqn. 9), and V_{water} is the volume of deionized water (Eqn. 10).

$$V_{stock\ collagen} = \frac{(V_{total} \times 1.2\ mg/mL)}{C_{stock\ collagen}} \quad \text{Eqn. 7}$$

$$V_{PBS} = \frac{V_{total}}{10} \quad \text{Eqn. 8}$$

$$V_{NaOH} = V_{stock\ collagen} \times 0.025 \quad \text{Eqn. 9}$$

$$V_{water} = V_{total} - (V_{stock\ collagen} + V_{PBS} + V_{NaOH}) \quad \text{Eqn. 10}$$

Note. The exact concentration of the collagen stock is lot dependent.

6. Transfer $V_{Cell\ stock}$ determined with Eqn. 2 to a 1.5 mL microcentrifuge tube and centrifuge at 1000 xg for 5 min.
7. Prepare the collagen suspension using the volumes calculated in Eqns. 7-10. First, combine the DI water, 10X PBS solution, and 1 M NaOH solution to an appropriate volume tube, mix by trituration, and then add the concentrated suspension of collagen I. Keep the prepared collagen solution on ice.

Note. The collagen suspension is viscous, pipette slowly to ensure transfer of the appropriate volume.

8. Mix the newly prepared collagen suspension thoroughly by trituration. Place a small volume (< 5 μ L) on a piece of pH paper. If the pH is not approximately 7.2, add 1 M NaOH until the correct pH is reached.
9. Aspirate the supernatant from the cell pellet.
10. Add $V_{Cell\ dep}$ calculated in Eqn. 3 to the cell pellet and mix until the suspension is homogenous.

Notes. Mixing is a critical step to ensure the same number of cells are deposited into each paper scaffold. Trituration, scraping the bottom and sides of the microcentrifuge tube, and mixing an XYZ-circular motion are all necessary to make a homogenous mixture. Take care not to introduce bubbles as they could harm the cells. Mix gently as not to kill the cells.

Take care to hold the tube by the top while mixing, as your hands will warm the collagen and accelerate gelation.

The mixing step is one of the first places to troubleshoot the deposition protocol if there is a large deviation between replicates.

11. Place 10 mL of fresh, 37 °C, 1X PBS into a 10 cm untreated petri dish.
12. Deposit the collagen-cell suspension into the paper scaffolds. For the small single zones, deposit 0.5 μ L with the 2.5 μ L micropipette. Place the pipette tip just above the center of the paper scaffold, slowly dispense the contents and allow it to wick into the paper and continue dispensing until you reach the first stop of the micropipette. Spread the collagen, moving in a circular motion from the center toward the wax borders. Ensure there are no 'dry' spots that contain no collagen. Once the volume has been spread evenly across the scaffold, dispense to the second stop of the pipette.

Note. When depositing the cell-collagen suspension, orient the paper perpendicular to the bench surface and hold it in place using a single "third hand" device, as demonstrated [Figure 1](#).

13. Allow the scaffold to dry for approximately 30 sec before transferring it to the Petri dish prepared in Step 12.

Note. Selecting an appropriate and consistent drying time is extremely important for obtaining reproducible datasets. If the scaffolds gelate for too long, they can dry out and reduce the overall cell viability. If the scaffolds do not gelate enough, the collagen might not fully crosslink, causing the beads/cells to fall out of the scaffolds.

The collagen is sufficiently dry when it is no longer shiny in appearance, but the cell-containing region is clearly still wet. Watch the first few scaffolds closely and choose a consistent drying time based on your observations. The small scaffolds require approximately 30 sec - 1 min to dry.

14. Repeat Steps 12 and 13 until 20 scaffolds have been deposited with cells.
15. Deposit three zones with 0.5 μ L of cell-free collagen. These scaffolds will serve as blanks to compare seeded fluorescence intensity and account for ECM autofluorescence.
16. Using the inverted microscope, check each scaffold to ensure the cells remain in the scaffold before continuing.

Note. The deposited cells look like clear spheres under the microscope, they should be packed throughout the scaffolds "culture area." Also, look at the medium for any cells that may have fallen out of the scaffolds. Some cells in the media is normal. If there are equal numbers of cells in the medium as there are in the paper scaffolds, then the gelation time was too short.

17. For the monolayer culture control, calculate the appropriate volume of medium (V_{medium}) and volume of cell stock ($V_{2D\ cell\ stock}$) required to prepare enough solution to deposit cell-containing medium along 20 wells of a 96-well plate with Eqns. 11 and 12. Each well will contain 4.0×10^4 cells and 200 μL (2.0×10^5 cells/ mL) cells of medium. Make ~ 5 mL of cell-containing media. C_{stock} is the same measured step 2.

$$V_{2D\ cell\ stock} = \frac{(2.0 \times 10^5 \frac{cells}{mL}) \times 5\ mL}{(C_{stock} \frac{cells}{mL})} \quad \text{Eqn. 11}$$

$$V_{medium} = 5\ mL - V_{2D\ cell\ stock} \quad \text{Eqn. 12}$$

18. In a 15 mL conical tube, add V_{medium} and $V_{2D\ cell\ stock}$. Inverter twice to mix cell containing solution. Pipette 200 μL of the cell containing into 20 wells of a fresh 96-well clear bottom plate. There should be 4.0×10^4 cells per well
19. Incubate the well plate and scaffolds overnight at standard culture conditions.

Measuring the fluorescence of the paper scaffolds

1. Transfer the cell-containing scaffolds to a fresh clear bottom 96-well plate with 100 μL of 1X DPBS. Using a pair of tweezers, carefully place the scaffolds on the surface of the DPBS, so they float at the air-liquid interface.

Note. The DPBS is used to ensure the scaffolds stay hydrated while preparing other scaffolds and the monolayer cells.

2. Aspirate the medium from the monolayer 96-well plate and the DPBS 3D 96-well plate just before analysis. If fluorescence measurements need to be taken sequentially (i.e., microscope, well plate reader), aspirate medium/DPBS from the plates one at a time.

Note. When aspirating from the paper scaffold-containing wells, do not disturb the paper scaffold floating at the air-liquid interface. A pipette can be used to remove DPBS slowly and gently

3. Measure the fluorescence intensity of each scaffold- and monolayer-containing well, and blank scaffolds or wells, with the flatbed scanner, well plate reader, or a fluorescence microscope. Collect a white light and fluorescent image of each well.

Note. For microscope capture, ensure the exposure and focus is the same for each scaffold.

Quantifying the number of viable cells in paper scaffolds using the CellTiter-Glo 2.0 assay

1. Aspirate any residual culture medium or DPBS from each well after measuring fluorescence and replace with 100 μ L of 1X DPBS. Add 100 μ L of 1X DPBS to blank scaffolds and wells.
2. Add 100 μ L of the CTG 2.0 reagent to each cell-containing scaffold and well. Add 100 μ L of CTG2.0 to blank scaffolds and wells.
3. Place the well plate on an XY orbital shaker and shake at room temperature for 15 min at about 750RPM.

Note. Ensure the medium is thoroughly mixed. Observe the scaffolds and well plates with an inverted microscope to ensure all the cells are lysed.

4. Transfer 100 μ L of the cell lysate solution to an opaque white 96-well plate. Analyze on a well plate reader with a luminescence readout.

Notes. If bubbles are present, use a 22G needle to pop them prior to analyzing the plate in the plate reader. To best visualize bubbles, Tilt plate 30° toward you and inspect along the top edge of each well.

We measure luminescence by capturing all photons in the range of 300-850 nm.

Analyzing the fluorescence of the paper scaffolds

1. Export images and analyze with desired analysis software (e.g., Fiji or ImageJ).
2. Inspect each scaffold to ensure there is a uniform distribution of fluorescence signal, and there are no regions of high and low (or no) fluorescence intensity, indicating a non-uniform spread during the deposition process. Visible differences in fluorescence are not acceptable, and the procedure (Protocols S2.2 a, b, and c) must be repeated.
3. Quantify the distribution of cells by obtaining a 3x3 matrix of fluorescence intensity on each scaffold. Calculate the average, standard deviation, and percent relative standard deviation (%RSD) of each scaffold. The %RSD within each scaffold should be less than 15%.
4. Quantify the average fluorescence intensity of the entire culture region in each scaffold to ensure reproducible deposition of cells. On a scaffold, make a circle object that is just smaller than the size of the cell culture zone. Using this circle, obtain the average intensity for all 20 scaffolds. Calculate the average, standard deviation, and %RSD of the dataset. The %RSD should be below 20%.

Note. Use the same-sized circle when obtaining intensity measurements between scaffolds.

5. To ensure the fluorescence intensity measurements are from the cells, calculate the signal-to-noise ratio. On a scaffold, make a circle object that is just smaller than the size of the cell culture zone, then calculate the average intensity of the cell-containing scaffolds and blank scaffolds. Divide the average fluorescence intensity of the cell-containing scaffolds by the average fluorescence intensity of the blank scaffolds. This value should be greater than 3.
6. Repeat this protocol on various days. The training set is complete when less than 20% %RDS is achieved on three separate depositing trials.

Analyzing the luminescence data to assess relative signal intensity from cells in the paper scaffolds using the CellTiter-Glo 2.0 assay

1. Export the luminescence data collected in above.
- 2a. Subtract the average monolayer blank luminescence signal from each of the cell monolayer samples.
- 2b. Subtract the average 3D blank luminescence signal from each of the 3D cell scaffold samples.
3. Calculate the average (\bar{I}), standard deviation, and %RSD for the background subtracted monolayer and 3D cell measurements.
4. Calculate the relative signal intensity of the 3D scaffolds (I_{Rel}), by dividing the average intensity of the 3D scaffolds (\bar{I}_{3D}) by the average intensity of the 2D scaffolds (\bar{I}_{2D}) multiplied by 100 percent, Eqn. 13.

$$I_{Rel} = \frac{\bar{I}_{3D}}{\bar{I}_{2D}} * 100\%$$

Eqn. 13

5. Evaluate the relative viability and the %RSD of the 3D scaffolds. Repeat this process until a relative viability of $\geq 85\%$ and a %RSD of $\leq 20\%$ is achieved on three separate depositing trials.